de	Gebrauchsanweisung	2
en	User's Manual	8
fr	Mode d'emploi	14
it	Istruzioni per l'uso	20
es	Instrucciones de empleo	26
pt	Manual de operação	32
ru	Руководство по применению	38

KaWe EUROLIGHT® E25 / E35 LED

CE

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User's Manual Ophthalmoscope KaWe EUROLIGHT® E25 / E35 LED, 5 apertures

Dear Customer, thank you for choosing a KaWe product. Our products are known for their high quality and long life. This KaWe product meets EC standards 93/42/EWG (standards for medical products).

Please read these instructions thoroughly and carefully before attempting to use this product and heed the given instructions for the care of the product!

Familiarize yourself fully with the operation of this product before attempting to use it.

Use: These ophthalmoscopes are only to be used by authorized professional personnel. They are only to be used for the purpose of examining the fundus of the eye.

Purpose: KaWe ophthalmoscopes are used to examine the fundus of the eye using direct ophthalmoscopy. The detection of characteristic changes in the fundus of the eye enables the diagnosis of numerous diseases. The conceptual design of the KaWe ophthalmoscope allows for its use without pupil dilating medication. In this case, however, the patient's field of sight is limited.

Unsuited use/contraindication of the products:

Any use other than that described here is not in accordance with the intended use of the unit. The manufacturer is not liable for any resulting damages. The user alone bears the risk. To power the KaWe EUROLIGHT®

E25, use exclusively the KaWe MedCenter® 5000 wall charging station.

KaWe EUROLIGHT® E25 / E35 LED, 5 apertures

-  **Slit aperture** for the detection of height fluctuations (from tumors or papilledema etc.)
-  **Small spot** for normal fundus examinations
-  **Large spot** for normal fundus examinations
-  **Red-free (green)- large spot RF** for detecting minute vein abnormalities, filters red light waves from the observation field
-  **Fixing agent** for detecting central and excentric fixation

Operation of the Ophthalmoscope

1. Reboss Lens: If both the observer and the patient are completely emmetropic (have perfect vision), it should be possible to obtain a clear fundus image without the use of a corrective lens. However, in most cases, both the observer's and the patient's eyes have refractive errors and in order to obtain a clear fundus image, it is necessary to correct the combination of these refractive errors. Hyperopic vision is corrected with a "plus" lens while myopic vision is corrected with a "minus" lens.

For example, consider a case where the observer's eye is myopic and the patient's eye is hyperopic.

(-5D): Correction of Observer's Myopia

(+3D): Correction of Patient's Hyperopia

→ (-2D): -4D Correction of Accommodation

In this case, the observer should use the -4D Rekoss lens to examine the patient's eyes in order to obtain a sharp fundus image. The power of the Rekoss lens is determined by the extent of the observer's and the patient's refractive errors and the distance between the ophthalmoscope and the patient's eye, etc.

The actual focussing is performed by rotating the Rekoss Disc while viewing the patient's fundus until

the observer achieves the sharpest focus of the fundus images. In case of high myopia or hyperopia, either a minus or plus auxiliary lens should be added.

The Diopter Indicator shows the combined dioptric power of the Rekoss and the auxiliary lenses. The red figure represents minus power, while the green figure represents plus power. The ranges of powers that are covered are shown below.

Position of the Auxiliary Lens Selectors	Diopter strength
•	-12D to +11D
+	+12D to +35D
-	-13D to -36D

2. Filter Lever: (a) Moving the Filter Lever to “F” inserts a 4000°K filter into the illumination path. This is used when illumination that is more like natural light is desired while observing the fundus of the eye. (b) When the Filter Lever is at the middle position, no filter is inserted into the illumination path. (c) Moving the Filter Lever to “P” inserts a polarizing filter into the illumination path.

3. Viewing Aperture Polarizing Filter: Viewing Aperture Polarizing Filter: Rotate the polarizing filter until the white mark is pointing upward, then set the Filter Lever to “P” and then while looking through the lens, ensure that the polarization axis is perpendicular to the polarization of the light. This serves to minimize reflection from the patient’s cornea. In order to obtain optimum balance between the brightness of the visual field and the amount of corneal reflection, the Viewing Aperture Polarizing Filter may be rotated. When the polarizing filter is in the locked position, the light coming into the viewing aperture is not polarized.

4. Aperture Dial: The Revolving Aperture Dial (7) inserts the following apertures into the illumination path:

(a) **Large spot:** This is used for ordinary fundus examinations.

(b) **Small spot:** This is used when viewing the fundus through a small pupil. It is especially useful when examining macula since the pupil tends to contract in response to bright light from the ophthalmoscope.

(c) **Red-free (green)- large spot:** Viewing the fundus with a red free filter makes the retinal blood vessels appear black against the grey background of

the fundus. This is particularly valuable for detecting minute superficial haemorrhages.

(d) **Fixing agent:** This aperture is used for the diagnosis of eccentric fixation. Direct a patient to look at the centre of the concentric circles while viewing the fundus. If an eccentric fixation is present, the patient’s fovea will not fall on the centre of the concentric circles.

(e) **Slit Aperture:** This aperture is used for the estimation of the height of various areas of the retina.

KaWe EUROLIGHT® E25

5. Assembling the KaWe MedCenter® 5000: To assemble the KaWe MedCenter® 5000 handles, push the head into the handle. Ensure that the center hole on the connector meets up with the fastening screw. Gently tighten the fastening screw. To remove, loosen the fastening screw and carefully pull out the head.

6. Operating the KaWe EUROLIGHT® E35 LED

(with click-connector): The ophthalmoscope head is designed to be clicked onto the battery handle. A locking device ensures a secure fit. To remove the head, push down the knurled sleeve. The lock will be released and the diagnostic head can be pulled off. Turn on by pressing and turning the rotary switch to the left. Turn the in-

strument off by turning the switch back to the right to its original position.

7. Changing the bulb: To replace the lamp, remove the head from the handle and pull out the old lamp. If necessary, clean the glass bulb of the new lamp with alcohol. The bulb should be free of fingerprints and grease. Insert the new lamp as follows: Align the side centring pin on the new lamp with the groove on the instrument head and slide the new lamp all the way into the socket.

Note: Proper operation of the instrument is only guaranteed when original KaWe vacuum bulbs are used. These bulbs ensure optimal illumination quality!

When handling the instrument, the examiner should try to keep the instrument as close as possible to his/her eye. If the examiner directs the light laterally into the fundal region at a distance from the subject of no more than 2-3 cm, a diagnosis can be made without the use of medication.

8. Window cover: When the instrument is not in use, close the cover on the viewing window to prevent dust from entering the instrument head.

How to View a Patient's Fundus: **1.** When examining the patient's right eye, hold the ophthalmoscope with your right hand and observe with your right eye. For examination of the left eye, use your left hand and eye. **2.** Put your index finger on the edge of the Rekoss Disc so that you can rotate it any time and hold the ophthalmoscope with your remaining four fingers. **3.** Place the ophthalmoscope on the side of your nose and position it so that you can look through the view-

ing aperture. **4.** Look through the viewing aperture and rotate the Rekoss disc so that you can clearly see an object placed at a distance. By doing this, your refractive error will be corrected. **5.** Seat the patient with his face toward the dark side of the room and direct him to look at an object placed at a distance. This procedure is necessary for keeping the patient's eye still and relaxed. Keeping the patient's eye still is essential to ophthalmoscopy since it is almost impossible to obtain a fundus image when the patient's eye is constantly moving. **6.** When examining the right eye, stand at the patient's right side. Hold the ophthalmoscope close to your right eye and direct the light beam into the patient's pupil at a distance of about 15 cm from the patient's eye. You will then be able to see a red reflection in the patient's pupil. **7.** Move slowly toward this red reflection while holding the ophthalmoscope close to your eye. Always keep the light beam in the patient's pupil while approaching the patient. When the ophthalmoscope is close enough to the patient's eye (2 to 3 cm), you will see the fundus image. **8.** If the fundus image is not in sharp focus, rotate the Rekoss disc and bring the image into focus. If it is not possible to get an image that is in focus within the range of powers of the Rekoss lenses, add a minus or a plus auxiliary lens by moving the Auxiliary Lens Selector and rotate the Rekoss disc again. **9.** When examining various areas of the fundus, hold the ophthalmoscope close to your eye and move it with your head around the patient's pupil, always keeping the same distance from the pupil. The light beam should always stay in the patient's pupil while moving. To stabilize the ophthalmoscope, put your middle finger on the patient's cheek and it will become easier to move the ophthalmoscope. **10.** If the reflection from the patient's cornea disturbs the fundus image, come as near to the patient as possible and direct the light beam through the lower

portion of the pupil while viewing the fundus through the upper portion of the pupil. If the corneal reflection still disturbs the observation, move the Filter Lever to "P". Then, rotate the Viewing Aperture Polarizing Filter so that its white mark is at the top. This position minimizes the corneal reflection. If you want to obtain a brighter visual field, you can do so by rotating the Viewing Aperture Polarizing Filter, but this will in turn increase the corneal reflection.

Further instructions, maintenance, storage:

When used and stored properly, this instrument will serve you reliably for many years. KaWe ophthalmoscopes are only to be used with KaWe-brand replacement parts and accessories.

Guarantee: When used under normal circumstances and with attention to these instructions (the User's Manual) we guarantee this product for two years after the date of purchase. Should you need further information or should your instrument require repair, please contact your dealer!

Warning: Never immerse the ophthalmoscope in any kind of liquid and always prevent all liquids from penetrating into the housing! Do not allow the ophthalmoscope to have direct contact with the eye. Long operation times can cause the parts of the instrument near the light to become very hot. The ophthalmoscope is not intended for use during operations and therefore simple cleaning and/or disinfection with an alcohol-based surface disinfecting agent is sufficient. Do not use abrasive cleaning agents! The ophthalmoscope cannot be sterilized.

KaWe ophthalmoscopes are to be used at room temperatures of not greater than 35° C. KaWe oph-

thalmoscopes are subject to special safety measures with regard to their electromagnetic compatibility. The operation of the KaWe ophthalmoscope can be influenced by portable and mobile high-frequency communication devices.

Workstation: Remove surface impurities with a paper towel or disposable cloth.

Storage and transport: No special requirements.

Cleaning preparation: No special requirements.

Manual cleaning: In order to avoid contact with dirt and dust, always store the ophthalmoscope in its packaging. A damp, soft and lint-free cloth can be used to clean the exterior of the instrument. The glass of the viewing window can be cleaned when necessary with a cotton swab moistened with a small amount of alcohol. Please avoid applying pressure on the glass while cleaning!

Automatic cleaning: Automated cleaning is neither necessary nor possible.

Disinfection: To disinfect the ophthalmoscope, wipe it with a cloth that has been moistened with alcohol.

Maintenance: Remove and replace damaged parts. If the instrument does not function properly, proceed as follows: **1.** Replace the light bulb, test functionality of instrument. **2.** Send the instrument in for repair.

Inspection and operation test: Check to ensure the smooth operation of all movable parts such as the aperture selection and correction lens dials. Turn on the instrument. Ensure that the head is properly

connected to the handle. For all instruments: visually check for damage and wear. Keep all electrical contact parts clean.

Packaging KaWe EUROLIGHT® E25: Packaged individually in standard polyethylene bags.

Packaging KaWe EUROLIGHT® E35 LED: Ophthalmoscope head and handle packaged in a hard case.

Sterilization: Not a surgical instrument, no sterilization necessary or possible.

Additional information: No special requirements

Storage: No special requirements.

Common accessories: Further information about these articles can be found on our homepage: **www.kawemed.de**

Manufacturer: **KaWe**

Contacting the manufacturer: Address or telephone number of the dealer or dial +49-7141-68188-0.

Symbol key:

	Manufacturer
	Date of manufacture
	Heed the User's Manual
	Separate disposal of electric and electronic devices
	Dispose of expired batteries separately at an appropriate collection center

	Complies with relevant EU guidelines
	Warning: electromagnetic field
	Only use device in dry rooms
	GOST-R certification for export articles to Russia
	Keep dry
	Temperature limits

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KIRCHNER & WILHELM GmbH + Co. KG