- de** Gebrauchsanweisung
- en** User's Manual
- fr** Mode d'emploi
- it** Istruzioni per l'uso
- es** Instrucciones de empleo
- pt** Manual de operação
- ru** Руководство по применению

Reflexhammer

BUCK

VARIOFLEX

Nadel
Needle
Piquenerf
Ago
Aguja
Agulha
Игла

Pinsel
Brush
Pinceau
Pennellino
Pincel
Pincel
Кисточка

de	Reflexhammer	4
en	Reflex hammer.....	8
fr	Marteau à réflexes	12
it	Martelletto per riflessi	16
es	Martillo de reflejo	20
pt	Martelo de reflexos	24
ru	Молоток неврологический	28

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Dear buyer! Thank you for choosing a KaWe product. Our products are known for high quality and durability. This KaWe device complies with the provisions of Regulation EU 2017/745 (European Regulation on Medical Devices), whereby it is classified as a medical product of class 1.

Please carefully consider this instruction and read all the provisions prior to use. Follow the maintenance recommendations.

Application: Only authorized qualified personnel can use tendon hammers and plexors.

Group of patients with indications for use: All

Intended use: Diagnostic tool for percussion (medical tapping).

Improper use of the device/contraindications for use: Any other use of use outside the above intended use shall be considered improper. The manufacturer is not responsible for damage resulting from such use. In this case, only the operator assumes all risks.

Maintenance of the tendon hammer and plexor: Tendon hammer is used to tap corresponding tendons freely, without effort in the wrist. This causes a reflex contraction of the corresponding muscles. Using a needle and a brush, one can check different skin reflexes (checking skin/surface sensitivity).

Additional information, maintenance and storage: The device will serve for years with proper handling and storage.

Secondary treatment limitations: Secondary treatment has a minor effect on these devices. The service life of the product is usually defined based on the degree of wear and use-related damage. Prompt secondary processing of the tool is recommended after its application.

Instructions

Transportation and storage: No special requirements.

Preparation for cleaning: No special requirements.

Manual cleaning: Wipe the tendon hammer with a damp cloth or a surface cleaner. Observe the concentrations specified by the manufacturer.

Disinfection: Disinfecting and cleaning solutions should be used in the correct dosage, according to the manufacturer. If the dosage is not enough, the solution is ineffective; too strong dosage causes damage to the tool. Observe the exposure time. Make fresh disinfecting and cleaning solutions regularly, as they lose their effectiveness.

Disinfection is the first processing procedure for contaminated items and materials in order to reduce the number of microorganisms and facilitate further cleaning. Avoid drying of dirt on the material. In addition, disinfection is aimed to protect personnel and the environment during the procedure. If possible, disinfect the tool at the place of its use as quickly as possible in one action.

Sterilization: Not applied for hammers! Only needles can be sterilized. The needles can be successfully sterilized with the following parameters: Fractional pre-vacuum (three times), sterilization temperature: min. 132 °C, max. 137 °C, exposure time: at least 3 min. (full cycle), drying time: 10 min.

Maintenance: Remove or replace the damaged parts.

Control and functional check: Check the ease of movement of moving parts, such as the needle and brush. Visually check for damages and wear. Before using, check the thread tightness and tight hold of elastic bands.

Package: Separate. A standard plastic bag can be used. The bag must be of sufficient size for the tool so that the clasp is not under tension.

Storage: Keep the tool away from dust, humidity and contamination during storage and transportation

Ambient temperature:

Operation: Room temperature

Transportation and storage: from $-10\text{ }^{\circ}\text{C}$ to $+50\text{ }^{\circ}\text{C}$

Relative air humidity from 10 % to 75%

Disposal: The product can be disposed of together with household garbage. Please observe the local regulations.

General accessories: To find out more about the products, visit our homepage: www.kawemed.com.

Manufacturer: **KaWe**

Basic UDI-DI: 4030155KaWe07SW

Manufacturer's contact details: Dealer's address or phone number: +49-7141-68188-0.

Warranty: The warranty for the device totals two years from the date of purchase with proper handling and following the User's Manual. If you have any additional questions or need repair, please, contact your dealer.

Symbol key:

	Attention.		CE conformity mark
	Manufacturer		Article number
	Lot code		Medical device
	Please follow the User's Manual		GOST-R certification of export goods imported into Russia

Information for users and patients:

All serious incidents involving this device must be immediately reported to the manufacturer and the competent authorities of the respective Member State in which the user and/or patient is located.

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KIRCHNER & WILHELM GmbH + Co. KG