ri-former®

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Instructions
Diagnostic Instruments

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


Description:

1. Openings at the rearside of Diagnostic station / Extension module.
2. Sliding cover.
3. Rocker switch
4. Green control lamp.
5. Application part / Handle.
6. Switching ring on the handle.
7. Handle holders.
8. Optional clock.

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Table of Contents

1. Important information to be observed before operation
2. Safety information and electromagnetic compatibility
3. Intended use
4. Attachment
5. Operation and function
6. Switching on and off
7. Cleaning and disinfection
8. Technical data
9. Electromagnetic compatibility

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1. Important information read prior to start-up

You have purchased a high quality **Riester** wall instrument, which has been manufactured according to the Directive 93/42 EEC and is subject to the strictest quality controls at all times. Read these instructions for use carefully before putting the unit into operation and keep them in a safe place. If you should have any questions, we are available to answer queries at all times. Our address can be found in these instructions for use. The address of our sales partner will be given upon request. Please note that all instruments described in these instructions for use are only to be used by suitably trained personnel. The perfect and safe functioning of this instrument is only guaranteed when original parts and accessories from **Riester** are used.

2. Safety information and electromagnetic compatibility

Meaning of the symbol on the model identification plate Attention, read accompanying papers!

Device of protective class II

Application part type B

Function earth

Protective grounding connection

MEDICAL – GENERAL MEDICAL EQUIPMENT
AS TO ELECTRICAL SHOCK, FIRE AND
MECHANICAL HAZARDS ONLY IN ACCORDANCE
WITH ANSI/AAMI ES60601-1 (2005) AND CAN/
CSA-C22.2 No. 60601-1 (2008)
Control No. 3TCH

Direct current

Alternating current

ON

OFF

Follow the instructions in the operation manual

Rel. Humidity

Used electrical and electronic devices should not be disposed of with normal household waste, but should be disposed of separately according to national or EU guidelines.

The instrument satisfies the requirements for electromagnetic compatibility. Please note that under the influence of unfavourable field strengths, e.g. during the operation of wireless telephones or radiological instruments, adverse effects on function cannot be excluded.

Attention!

There is a possible danger of inflammation of gases, if the instrument is operated in the presence of inflammatory mixtures or mixtures of pharmaceuticals and air or oxygen or laughing gas! Never attempt to take the instrument apart! There is a danger of life-threatening electrical shock. Unplug the instrument before cleaning or when disinfecting.

For devices of protection class I : „Attention : To avoid the risk of an electric shock the device must only be connected to the public power supply with a protective earthing conductor“.

3. Intended use

The Wall /Mobile, Anesthesia device **ri-former®** described in these instructions for use was manufactured for use with various instrument heads and modular components for non-invasive diagnostics.

4. Attachment

a.) Drilling instructions/drilling plan

The drilling instructions and the drilling plan are enclosed separately. Follow the drilling instructions in order to drill the holes in the wall.

b.) Attaching the wall mounting plates

After you have drilled the holes, take the plugs supplied and push them into the holes as far they will go. Take the wall mounting plate and hold it onto the wall so that the screws can be pushed through the holes of the mounting plate into the plugs. Now screw in the screws with a screw driver, as far as they will go.

c.) Attachment of the diagnostic station

When all screws have been screwed in tightly, take the diagnostic station and guide the screw heads through the openings (1). Then press the diagnostic station downwards until it snaps into place.

d.) Attachment of the extension module

Connect the diagnostic station and the extension module with **the** help of the connecting cable. In order to plug in the connecting cable, remove the sliding cover (2) of the diagnostic station. Close the casing opening of the extension module, which is not needed, with the sliding cover (2). Take the extension module and guide the screw heads through the openings (1). Then press the extension module downwards.

Attention:

Take care that the connecting cable does not get caught behind the extension module. Push the connecting cable into the groove provided on the reverse side of the extension module.

4.1 ri-former® mobile

Please follow the assembly instructions for the mobile stand with **ri-former®**. The assembly instructions is included in the mobile stand delivery package.

4.2 ri-former® Anesthesia**Assembly of the universal clamp.**

Please check if the designated wall rail is rigidly mounted on the wall. Fix the universal clamp on the determined place of the wall rail and tighten very well the locking screw. Put the pre-assembled **ri-former®** anesthesia device on the universal clamp and insert it. Please ensure that both pins are introduced in the universal clamp. Afterwards tighten the **ri-former®** anesthesia device with the lateral screw.

5. Operation and function**Putting the diagnostic station into service with or without extension module:**

Put the plug into the electrical socket. The optional clock (8) starts to blink. You can adjust it to local time by repeatedly pressing the keys; with the left key marked HR and the right key marked MIN. Move the handle (5) upwards out of the handle holder (7) and attach the desired instrument head by placing it with the two projecting guide cams onto the handle. Press the instrument head lightly onto the handle and turn the handle in a clockwise direction until it stops. Removal of the instrument head is carried out by turning in a counter-clockwise direction.

6. Switching on and off

Switch on the instrument by using the rocker switch (3). The green control lamp (4) in the rocker switch (3) indicates that the instrument is ready to use. Each handle (5) is automatically ready to operate at 100% light intensity as soon as it is taken out of the handle holders (7). The handle is switched off automatically by putting back into the handle holder. The handle (5) is automatically switched off when replaced back into the handle holder (7).

6.1 rheotronic® for light intensity modulation

the modulation of light intensity can be done with the handle; you only have to tip the switching ring clockwise direction or against clockwise direction and the light gets stronger or weaker.

Attention!

The handle gets off automatically after abt. 3 minutes. Make sure that no more than 3 handles (5) are used at the same time! If more than 3 handles are used at the same time, the transformer in the instrument may become overloaded and switch itself off.

7. Care instructions:**General Note**

Cleaning and disinfection of medical products are meant to protect patients, users, and third persons and lead to value retention of medical products. Due to product design and materials used, there is no possibility to define the maximum limit of re-processing cycles. Lifetime of medical products is determined by its function and gentle use. Before sending back defective products for repair, the described re-processing cycles have to be applied

and followed.

Cleaning and disinfection

The diagnostic station **ri-former**[®] with extension module can be cleaned on the outside with the exception of the clock glass cover (8), using a damp cloth until optical cleanliness is achieved. Wipe disinfection only according the instructions of use of respective manufacturer of disinfectant. Only disinfectants with proven effectiveness in consideration of national guidelines should be used. After disinfection, please wipe the instruments using a damp cloth in order to eliminate eventual remains of the disinfectants.

Attention!

- We recommend unplugging the instrument before cleaning or disinfection.
- Take care while cleaning and disinfecting that no liquid enters inside the instrument.
- Place removable parts of the diagnostic station **ri-former**[®] with expansion modules (handle, instrument heads) never in liquids!
- The devices are not meant to undergo machine-processed maintenance and sterilization. It may lead to irretrievable damages!

8. Technical data

Model:	Voltage supply Diagnostic station ri-former [®]
Voltage:	120V~/50-60 Hz 230V~/50-60 Hz 240V~/50-60 Hz
Input:	31 VA
Output:	1 x 3,5Vdc / 2 x 12Vac
Output current:	1 x 700mA / 2 x 1000mA
Classification:	Application part type B
Working temperature:	0° C to + 40° C, 10% up to 85 % relative humidity
Storage location:	-5° C to + 50° C, 10% up to 85 % relative humidity
Airpressure:	700 bis 1050hPa
Dimensions	
Diagnostic station:	200 x 180.5 x 75 mm
Extension module:	200 x 100 x 75 mm
Weight Diagnostic station:	1450 g
WeightExtension module:	490 g
Switn-on time:	ON: 1 Min OFF: 5 Min

9. Electromagnetic compatibility

Medical electrical equipment is subject to special precautionary measures with regard to electromagnetic compatibility (EMC). Portable and mobile high-frequency communication equipment can influence medical electrical equipment. This ME device is intended for operation in an electromagnetic environment as specified below. The user of the device should ensure that it is operated in such an environment. The ME device must not be used directly next to or arranged in a stack with other devices. If the device has to be operated near to or in a stacked arrangement with other devices, then the ME device should be monitored in order to verify that it operates as intended in this arrangement.

This ME device is intended exclusively for use by professional medical staff. This device can cause radio interference and can disrupt the operation of equipment nearby. Suitable remedial measures, such as for instance re-alignment, re-arrangement of the ME device or shielding, can become necessary.

Guidelines and manufacturer's declaration - electromagnetic emissions		
The Diagnostic station ri-former® is intended for operation in an electromagnetic environment as specified below. The customer or the user of the Diagnostic station ri-former® should ensure that it is used in such an environment.		
Emission measurements	Compliance	Electromagnetic environment - guidelines
HF emissions according to CISPR 11	Group 1	The Diagnostic station ri-former® employs HF energy solely for an internal function. Its HF emission is therefore very low and it is unlikely that neighboring electronic devices will be affected by interference.
HF emissions according to CISPR 11	Class B	The Diagnostic station ri-former® is intended for use in all facilities, including living quarters and such as are directly connected to a public power supply that also supplies buildings that are used for residential purposes.
Harmonics emissions according to EC61000-3-2	Not applicable	
Voltage fluctuation / flicker emissions according to IEC61000-3-3	Not applicable	

Guidelines and manufacturer's declaration - electromagnetic immunity

The Diagnostic station **ri-former**[®] is intended for operation in an electromagnetic environment as specified below. The customer or the user of the Diagnostic station **ri-former**[®] should ensure that it is used in such an environment.

Immunity tests	IEC 60601 test level	Compliance level	Electromagnetic environment -guidelines
Electrostatic discharge (ESD) according to IEC61000-4-2	± 6 kV contact discharge ± 8 kV air discharge	± 6 kV contact discharge ± 8 kV air discharge	Floors should be made of wood or concrete or be covered with ceramic tiles. If the floor is covered with a synthetic material, the relative air humidity must be at least 30%.
Fast transient electrical interference/ bursts according to IEC61000-4-4	± 2 kV for power lines ± 1 kV for input and output lines	± 2 kV for power lines ± 1 kV for input and output lines	The quality of the supply voltage should correspond to that of a typical business or hospital environment.
Surges IEC61000-4-5	± 1 kV voltage phase-to phase ± 2 kV voltage phase-to earth	± 1 kV voltage phase-to phase ± 2 kV voltage phase-to earth	The quality of the supply voltage should correspond to that of a typical business or hospital environment.
Voltage dips, short-time interruptions and fluctuations in the supply voltage according to IEC61000-4-11	$\leftarrow 5\% U_T$ ($\rightarrow 95\%$ drop in U_T) for 0.5 cycles $40\% U_T$ (60 % drop in U_T) for 5 cycles $70\% U_T$ (30 % drop in U_T) for 25 cycles $\leftarrow 5\% U_T$ ($\rightarrow 95\%$ drop in U_T) for 5 s	$\leftarrow 5\% U_T$ ($\rightarrow 95\%$ drop in U_T) for 0.5 cycles $40\% U_T$ (60 % drop in U_T) for 5 cycles $70\% U_T$ (30 % drop in U_T) for 25 cycles $\leftarrow 5\% U_T$ ($\rightarrow 95\%$ drop in U_T) for 5 s	The quality of the supply voltage should correspond to that of a typical business or hospital environment. If the user of the Diagnostic station ri-former [®] requires continued operation during power mains interruptions, it is recommended that the Diagnostic station ri-former [®] be powered from an uninterruptible power supply or a battery.
Magnetic field at the mains frequency (50Hz) according to IEC61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at level characteristic of a typical location in a typical commercial or hospital environment

Note - U_T is the alternating supply voltage prior to application of the test level.

Guidelines and manufacturer's declaration - electromagnetic immunity

The Diagnostic station **ri-former®** model is intended for operation in the electromagnetic environment specified below. The customer or the user of the Diagnostic station **ri-former®** should ensure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidelines
<p>Conducted HF interference according to IEC 61000-4-6</p> <p>Radiated HF interference according to IEC 61000-4-3</p>	<p>3 Vrms 150 kHz to 80MHz</p> <p>3 V/m 80 MHz to 2.5GHz</p>	<p>3 V effective-value</p> <p>3 V/m</p>	<p>Portable and mobile radio equipment should not be used within a distance from the Diagnostic station ri-former® including cables, that is less than the recommended safety distance as calculated by the equation that is appropriate for the transmission frequency.</p> <p>Recommended safety distance:</p> <p>$d = 1.2\sqrt{P}$</p> <p>$d = 1.2\sqrt{P}$ 80 MHz to 1000 MHz</p> <p>$d = 2.3\sqrt{P}$ 1400 MHz to 2.5 GHz</p> <p>where P is the nominal power of the transmitter in Watts (W) as specified by the manufacturer of the transmitter, and d is the recommended safety distance in meters (m). The field strength of stationary radio transmitters should be less than the compliance level b at all frequencies as verified by an on-site test a</p> <p>Interference is possible in the vicinity of equipment marked with the following symbol</p> 

Note 1: At 80 MHz and 800 MHz, the higher value applies.

Note 2: These guidelines may not apply in all situations. The propagation of electromagnetic waves is influenced by reflection and absorption by buildings, objects and people.

▪ The field strength of stationary transmitters, such as base stations of wireless telephones and mobile field radio services, amateur radio stations, AM and FM radio and television transmitters cannot be precisely determined theoretically in advance. In order to determine the electromagnetic environment due to stationary HF transmitters, an investigation of the location is advisable. If the field strength determined at the location of the Diagnostic station **ri-former**[®] exceeds the compliance level indicated above, then the Diagnostic station **ri-former**[®] must be monitored with regard to its normal operation at each place where it is used. If unusual performance characteristics are observed, additional measures such as re-alignment of the Diagnostic station **ri-former**[®] or its removal to another place may be necessary.

▪ In the frequency range of 150 kHz to 80 MHz, the field strength should be smaller than 3 V/m.

Recommended safety distances between portable and mobile HF communication devices and the Diagnostic station **ri-former**[®]

The Diagnostic station **ri-former**[®] is intended for operation in an electromagnetic environment in which the radiated HF interference is monitored. The customer or user of the Diagnostic station **ri-former**[®] can help prevent electromagnetic interference by observing minimum distances between portable and mobile HF communication equipment (transmitters) and the Diagnostic station **ri-former**[®] as recommended below, depending on the maximum output power of the communication equipment.

Nominal power of the Transmitter W	Safety distance that applies to the transmitter frequency m		
	150 KHz to 80 MHz d = 1.17VP	80 MHz to 800 MHz d = 1.17VP	800 MHz to 2.5 GHz d = 2.33VP
0.01	0,11667	0,11667	0.23333
0.1	0,36894	0,36894	0.73786
1	1,16667	1,16667	2.33333
10	3,68932	3,68932	7.37865
100	11.66667	11.66667	23,33333

For transmitters whose nominal power is not indicated in the table above, the distance can be determined using the equation belonging to the respective column, where P is the nominal power of the transmitter in Watts (W) as specified by the manufacturer of the transmitter.

Note 1: At 80 MHz and 800 MHz, the distance for the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. The propagation of electromagnetic waves is influenced by reflection and absorption by buildings, objects and people.

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MEDICAL — GENERAL MEDICAL EQUIPMENT AS TO ELECTRICAL SHOCK, FIRE AND MECHANICAL HAZARDS ONLY IN ACCORDANCE WITH ANSI/AAMI ES60601-1 (2005) AND CAN/CSA- C22.2 No. 60601-1 (2008) Control No. 3TCH

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ri-scope® L otoscope

1. Purpose

The **Riester** otoscope described in these operating instructions is produced for illumination and examination of the auditory canal in combination with **Riester** ear specula.

2. Fitting and removing ear specula

Either **Riester** disposable ear specula (blue colour) or reusable **Riester** ear specula (black colour) can be fitted to the otoscope head. The size of the ear specula is marked at the back of the speculum.

L1 and L2 otoscopes

Screw the speculum clockwise until noticeable resistance is felt. To remove the speculum, screw the speculum counter clockwise.

L3 otoscope

Fit the chosen speculum on the chrome-plated metal fixture of the otoscope until it locks into place. To remove the speculum, press the blue ejection button. The speculum is automatically ejected.

3. Swivel lens for magnification

The swivel lens is fixed to the device and can be swivelled 360°.

4. Insertion of external instruments into the ear

If you wish to insert external instruments into the ear (e.g. tweezers), you have to rotate the swivel lens (approx. 3-fold magnification) located on the otoscope head by 180°. Now you can use the operation lens.

5. Pneumatic test

To perform the pneumatic test (= examination of the eardrum) you require a ball, which is not included in the normal delivery package, but can be ordered separately. The tube for the ball is attached to the connector. Now you can carefully insert the necessary volume of air into the ear canal.

ri-scope® L ophthalmoscope

1. Purpose

The **Riester** ophthalmoscope described in these operating instructions is produced for the examination of the eye and the eyeground.

2. Lens wheel with correction lens

The correction lens can be adjusted on the lens wheel. The following correction lenses are available:

L1 and L2 ophthalmoscopes

Plus: 1-10, 12, 15, 20, 40.

Minus: 1-10, 15, 20, 25, 30, 35.

L3 ophthalmoscope

Plus: 1-45 in single steps

Minus: 1-44 in single steps

The values can be read off in the illuminated field of view. Plus values are displayed in green numbers, minus values with red numbers.

3. Apertures

The following apertures can be selected with the aperture hand-wheel:

L1 ophthalmoscope

Semi-circle, small/medium/large circular aperture, fixation star, slit and red-free filter.

L2 ophthalmoscope

Semi-circle, small/medium/large circular aperture, fixation star and slit.

L3 ophthalmoscope

Semi-circle, small/medium/large circular aperture, fixation star, slit and grid.

Aperture

function

	Small circle:	to reduce reflection for small pupils
	Medium circle:	
	Semi-circle:	
	Large circle:	for normal examination results
	Karo:	for topographic determination of retina changes
	Light slit:	to determine differences in level
	Fixation star:	to ascertain central or eccentric fixation

4 Filters

Using the filter wheel, the following filters can be switched for each aperture:

L1 ophthalmoscope

The L1 instrument head is supplied without a filter wheel. (the red filter is contained in the aperture wheel)

L2 ophthalmoscope

Red-free filter, blue filter and polarisation filter.

L3 ophthalmoscope

Red-free filter, blue filter and polarisation filter.

Filter

function

Red-free filter:	contrast enhancing to assess fine vascular changes, e.g. retinal bleeding
Polarisation filter:	for precise assessment of tissue colours and to avoid retinal reflections
Blue filter:	for improved recognition of vascular abnormalities or bleeding, for fluorescence ophthalmology

For L2 + L3, every filter can be switched to every aperture.

5. Focussing device (only with L3)

Fast fine adjustment of the examination area to be observed is achieved from various distances by turning the focussing wheel.

ri-scope® retinoscope (skiascope)

1. Intended use

The **ri-scope®** retinoscope Slit and **ri-scope®** retinoscope Spot described in these operating instructions (also called skiascopes) have been manufactured for examining the refraction of the eye (refractive error).

2. Function

Rotation and focusing of the slit and/or spot image may now be effected by the knurled screw.

3. Rotation

The slit or spot image may be rotated by 360° by the control. Each angle may be directly read from the scale on the retinoscope.

4. Fixation cards

Fixation cards are suspended and fixed on the object side of the retinoscope into the bracket for the dynamic skiascope.

5. Slit/Spot design

The slit retinoscope may be converted to a spot retinoscope by exchanging the slit lamp against a spot lamp.

ri-derma dermatoscope

1. Intended use

The **ri-derma** dermatoscope described in these operating instructions has been produced for early recognition of melanotic skin changes (malign melanoma).

2. Focusing

Focus the magnifying glass by rotating the eyepiece ring.

3. Skin adapters

Two skin adapters are supplied:

1) Including a scale of 0 - 10 mm for measuring melanotic skin changes, such as malign melanoma. = article number 10969

2) without scale = article number 10968

Both skin attachments can be removed easily and exchanged.

ri-scope® F.O. bent arm illuminator

1. Intended use

The bent arm illuminator described in these instructions for use was manufactured for illuminating the mouth and pharynx.

ri-scope® F.O. nasal speculum

1. Intended use

The nasal speculum described in these instructions for use was manufactured for illumination and examination of the inside of the nose.

2. Function

Two types of operation are possible:

a) Quick retraction

Press down the adjusting screw on the instrument head with

the thumb. In this adjustment, the position of the shank of the speculum cannot be changed.

b) Individual retraction

Turn the adjusting screw in a clockwise direction until you obtain the desired speculum opening. The shanks close again when the screw is turned in a counter-clockwise direction

3. Swivel lens

A swivel lens with an approx. 2.5-fold magnification is to be found on the nasal speculum, which can be simply removed or replaced again in the opening provided in the nasal speculum.

ri-scope® F.O. tongue blade holder

1. Intended use

The tongue blade holder described in these instructions for use was manufactured for the examination of the mouth and throat in combination with commercially available wooden and plastic blades.

2. Function

Insert a commercially available wooden or plastic tongue blade through the opening below the light outlet until it stops. After the examination, the blade can be easily removed by pushing the ejector.

ri-scope® human operations otoscope without speculum

1. Intended use

The **Riester** operation otoscope described in these instructions for use was manufactured for the illumination and examination of the auditory canal as well as for small operations in the auditory canal.

2. Attachment and removal of ear specula for human medicine

Place the desired speculum onto the black holder of the operation otoscope so that the recess on the speculum fits into the guide of the holder. Fix the speculum by turning it in a counter-clockwise direction.

3. Swivel lens for magnification

There is a small magnification lens which can be swivelled 360° on the operation otoscope with approx. 2.5-fold magnification.

4. Insertion of external instruments into the ear

The operation otoscope has been designed so that external instruments can be inserted into the ear without a problem.

ri-scope® Veterinary operation otoscope without speculum

1. Intended use

The **Riester** operation otoscope described in these instructions for use was manufactured solely for use in animals and thus for veterinary medicine. It can be used for illumination and examination of the auditory canal as well as for small operations in the auditory canal.

2. Attachment and removal of ear specula for veterinary medicine

Place the desired speculum onto the black holder of the operation otoscope so that the recess on the speculum fits into the guide of the holder. Fix the speculum by turning it in a counter-clockwise direction.

3. Swivel lens for magnification

There is a small magnification lens which can be swivelled by 360° on the operation otoscope with approx. 2.5-fold magnification.

13. Replacing the lamp

L1 otoscope

Remove the specula fitting from the otoscope. Screw out the lamp counter clockwise.

Screw in the new lamp clockwise and replace the specula fitting. L2, L3 otoscopes, ri-derma, bent-arm illuminator, nasal speculum and blade holder

Screw the instrument head off the battery holder. The lamp is located at the base of the instrument head. Pull the lamp out of the instrument head with thumb and forefinger or a suitable tool. Insert a new lamp.

Ophthalmoscopes

Remove the instrument head from the battery holder. The lamp is located at the base of the instrument head. Remove the lamp from the instrument head with thumb and forefinger or a suitable tool. Insert a new lamp.

Caution: The pin on the lamp must be inserted into the guide groove on the instrument head.

Veterinary/human operation otoscope

Screw the lamp out of the fixture in the operation otoscope and screw in a new lamp.

Instrument heads: Retinoscope slit and spot

Remove the instrument head from the battery handle. The lamp is located in a sleeve at the base of the instrument head. Remove the lamp from the sleeve using the thumb and index finger or a suitable tool. Insert the new lamp firmly into the sleeve and replace the sleeve back into the instrument head so that the base of the lamp fits into the slot on the instrument head.

Information on care, cleaning and disinfection

All **Riester ri-scope®** instrument heads can be cleaned on the outside with a moist cloth. Furthermore, the following disinfectants can be used: Aldehydes (formaldehyde, glutaraldehyde, aldehyde fission products) or tensides. All parts of instruments with the exception of the swivel lens, magnifying glass and the cover glass can also be disinfected with alcohols. When using these substances it is absolutely essential to follow the instructions of the manufacturer. A soft and, as far as possible, lint-free cloth or cotton bud can be used as an auxiliary aid for cleaning or disinfection.

Attention

Never place the instrument heads in liquid. Take care that no liquid enters inside the casing.

a) Sterilisation

According to the current school of thought (Test Centre for Medical Devices in Tübingen), sterilisation is only prescribed in the case of operative procedures.

b) Reusable ear specula


Although sterilisation is not necessary as described in a), it is nevertheless possible. The reusable ear specula can be sterilised at 134°C and 10 minutes retention time in a steam sterilizer.

Putting the instruments heads into operation

Place the desired instrument head onto the attachment on the handle so that the two recesses on the lower part of the instrument head sit on top of the two projecting guide cams of the battery handle. Press the instrument head lightly onto the handle and turn the handle in a clockwise direction until it stops. To remove the head turn it in a counter-clockwise direction.

Putting the anti-theft security into operation

1 Function

Place the desired instrument head onto the attachment on the handle so that the two recesses on the lower part of the instrument head sit on top of the two projecting guide cams of the battery handle. Press the instrument head lightly onto the handle and turn the handle in a clockwise direction until it stops.

In order to activate the anti-theft security, turn the Allen screw (b) using the Allen key (a) (included with the instrument head) until it stops. The instrument head can now no longer be removed from the handle. In order to deactivate the anti-theft security, the Allen screw (b) has to be unscrewed again using the Allen key (a).

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WARRANTY

This product has been manufactured under the strictest quality standards and has undergone a thorough final quality check before leaving our factory. We are therefore pleased to be able to provide a warranty of **2 years from the date of purchase** on all defects, which can verifiably be shown to be due to material or manufacturing faults. A warranty claim does not apply in the case of improper handling. All defective parts of the product will be replaced or repaired free of charge within the warranty period. This does not apply to wearing parts. For **R1 shock-proof**, we grant an additional warranty of 5 years for the calibration, which is required by CE-certification. A warranty claim can only be granted if this Warranty Card has been completed and stamped by the dealer and is enclosed with the product. Please remember that all warranty claims have to be made during the warranty period. We will, of course, be pleased to carry out checks or repairs after expiry of the warranty period at a charge. You are also welcome to request a provisional cost estimate from us free of charge. In case of a warranty claim or repair, please return the **Riester** product with the completed warranty Card to the following address:

**Rudolf Riester GmbH
Dept. Repairs RR
Bruckstr. 31
72471 Jungingen
Germany**

Serial number or batch number
Date, Stamp and signature of the specialist dealer

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